

Meditec 1020

Portable Emergency Ventilator



TECHNICAL SPECIFICATIONS

Pneumatically Driven Ventilator	✓
Driven from pipeline or cylinder operated with a non-return valve	✓
IPPV Ventilation Mode	✓
Adult & Pediatric	✓
I:E Ratio	Preset (1:2)
Oxygen Concentration	100%
Minute Volume with 100% O ₂	2–15 Ltr/min
Breathing Frequency (Adjustable)	8–30 bpm
Noise	<65 db
Cylinder Time - 3 Ltr. cylinder when fully charged with MV of 6 l/m)	>60 min
Driving Gas Pressure	280–600 kPa
Dimensions (mm) W x H x D	116 x 61 x 191
Weight	1.5 Kg
Storage Temperature	0°C–50°C
Operating Temperature	0°C–50°C
Storage Atmospheric Pressure	50–110kPa
Operating Atmospheric Pressure	70–110kPa

STANDARD ACCESSORIES

- Ventilator Main Unit
- High Pressure Hose
- Breathing Hose
- 2/3 Ltr. O₂ Cylinder
- Breathing Valve
- Face mask
- Regulator cum Flowmeter
- Cylinder Refilling Attachment
- Dust Proof Carry Bag



Dust Proof Carry Bag



Meditec International England Limited

8, Pinner View, Harrow, Middx,

HA1 4QA, United Kingdom

Tel: +44 203 6423808

E-mail : sales@meditecengland.co.uk

Website : www.meditecengland.co.uk

*In view of continued product improvement Meditec International England Limited reserves the right to make changes in the equipment which may affect the information contained in this leaflet without notice.